

Lensing Keys:

Reduce: [

Increase:]

Zoom In: .

Zoom Out: ,

No zoom: /

Ruler: ;

Default: \

Previous lens: ↑ [, {

Next lens: ↑] , }

None: ↑ \ , |

Put-down/pick-up: p

Snapshot: ↑ d, D

Lens Shape: ↑ l, L

Toggle visible/hidden: l

Toggle magnification modes: m

Histosearch lens: increase sensitivity with k/reduce with j, trigger search with h, delete current selection with x.