



LigPlot+ for Mac Users

1. Visit: <https://www.ebi.ac.uk/thornton-srv/software/LigPlus/>
2. In the menu on the left, click “Academic License”
3. Fill in all required fields. After clicking “Accept”, you should receive an email confirmation shortly
4. Click on the link in the email and fill in your email and given password
5. Ensure that you have Java Sun environment downloaded:
 - a. <https://www.java.com/en/download/manual.jsp> Click the “Mac x64” link
6. When you finish downloading LigPlus, find the zip folder in finder



a.

7. Double click on the zip file, and move the resulting LigPlus folder into the Users folder:



8. Once you know that LigPlus is in your Users folder, open up “Terminal”

- a. Click on the magnifying glass in the top bar next to the WiFi symbol, and search “Terminal”
9. Copy and paste the following into terminal:
 - a. `java -cp /Users/ -jar /Users/LigPlus/LigPlus.jar`
 - b. When you first open LigPlus, you will be brought to a window titled “Paths and Directories”
 - c. Here, ensure that in Temporary Directory, it says “/tmp” and NOT “C:\tmp”

After downloading the proper Miniconda installer from Github link, go to downloads folder and make a New Terminal at Folder:

Run: `bash Miniconda3-latest-MacOSX-x86_64.sh`