

## Overview of the voices available in AsTeRICS GRID

### iPad

- Recommended browser: Firefox
- Online and offline voices made available by Apple are listed in the drop-down list of voices in AsTeRICS GRID. The quality varies.
- A few of these voices are listed by Apple as “male” but don’t sound deep. It might be a mistake by Apple or possibly a male teen before voice break.
- The current system voice (Accessibility > Spoken Content > Voice) is available offline and is mostly high quality. It doesn’t get listed in the dropdown but activates when the voice is set to default. The iPad language must be set to that language (General > Language & Region) for the voice to work. If multiple voices for that language have been downloaded, the one first listed in the settings will be used. Changing the order in the settings is not possible, so only download the voice you want to use.
- Since only the current system voice is available, it’s not possible to combine system voices for multilingual users, but it’s possible to combine one system voice with the offline and online voices from the browser

Language	Female voice	Male Voice	Child Voice	system / offline / online	Name
Arabic	x	x		system	
Bangla	x			system	
Basque	x			system	
Bhojpuri	x			system	
Bulgarian	x			system	
Catalan	x	x		system	
Chinese, Hongkong (Cantonese)	x	x		system	
Chinese, Mandarin	x	x		system	
Chinese, Shanghainese	x			system	
Croatian	x			system	
Czech	x			system	
Danish	x	x		system	
Dutch	x	x		system	
English, Australia	x	x		system	
English, India	x	x		system	
English, Ireland	x	x (Siri)		system	
English, Scotland	x			system	
English, South Africa	x	x (Siri)		system	
English, UK	x	x		system	
English, US	x	x		system	

Finnish	x	x		system	
French	x	x		system	
Galician	x			system	
German	x	x		system	
Greek	x	x		system	
Hebrew	x	x (Siri)		system	
Hindi	x	x		system	
Hungarian	x			system	
Indonesian	x			system	
Italian	x	x		system	
Japanese	x	x		system	
Kannada	x			system	
Kazakh	x			system	
Korean	x	x		system	
Lithuanian	x			system	
Malay	x			system	
Marathi	x			system	
Norwegian, Bokmål	x	x		system	
Persian		x		system	
Polish	x	x		system	
Portuguese	x	x		system	
Portuguese, Brazil	x	x		system	
Romanian	x			system	
Russian	x	x		system	
Slovak	x			system	
Slovenian	x			system	
Spanish	x	x		system	
Spanish, Argentina	x	x		system	
Spanish, Chile	x			system	
Spanish, Colombia	x	x		system	
Spanish, Mexico	x	x		system	
Swedish	x	x		system	
Tamil	x			system	
Telugu	x			system	
Thai	x			system	
Turkish	x	x		system	
Ukrainian	x			system	
Vietnamese	x			system	
Arabic		x		offline	Majed
Bulgarian	x			offline	Daria
Catalan	x			offline	Montse
Chinese	x			offline	Tingting
Chinese, Hong Kong	x			offline	Sinji

Chinese, Taiwan	x			offline	Meijia
Croatian	x			offline	Lana
Czech	x			offline	Zuzana
Danish	x			offline	Sara
Dutch		x		offline	Xander
Dutch, Belgium	x			offline	Ellen
English, Ireland	x			offline	Moira
English, South Africa	x			offline	Tessa
Finnish	x			offline	Satu
French		x		offline	Thomas, Jacques
Greek	x			offline	Melina
Hebrew	x			offline	Carmit
Hindi	x			offline	Lekha
Hungarian	x			offline	Tünde
Indonesian	x			offline	Damayanti
Japanese	x			offline	Kyoko
Malay	x			offline	Amira
Norwegian, Bokmål	x			offline	Nora
Polish	x			offline	Zosia
Portuguese	x			offline	Joana
Romanian	x			offline	Ioana
Russian	x			offline	Milena
Slovak	x			offline	Laura
Swedish	x			offline	Alva
Thai	x			offline	Kanya
Ukrainian	x			offline	Lesya
Vietnamese	x			offline	Linh
Afrikaans			male teen?	online	
Albanian		x		online	
Arabic	x	x		online	
English, Australia	x	x		online	
Bangla, Bangladesh	listed as female, sounds male	x		online	
Bangla, India	x	x		online	
Bosnian		x		online	
Catalan			male teen?	online	
Chinese	x	x		online	
Chinese, Hong Kong	x			online	

Chinese, Hong Kong		x		online	
Chinese, Taiwan	x	x		online	
Croatian		x		online	
Czech	x			online	
Danish	x			online	
Dutch	x	x		online	
English, UK Fallback	x			online	
English, UK	x	x		online	
English, US	x	x		online	
Estonian		x		online	
Filipino	x			online	
Finnish	x			online	
French	x	x		online	
French Canadian	x	x		online	
German	x	x		online	
Greek	x			online	
Hindi	x	x		online	
Hungarian	x			online	
Icelandic		x		online	
Indonesian	x	x		online	
Italian	x	x		online	
Japanese	x	x		online	
Korean	x	x		online	
Latin		x		online	
Latvian		x		online	
Moldavian	x			online	
Montenegrin			male teen?	online	
Nepali	x			online	
Norwegian	x	x		online	
Polish	x	x		online	
Portuguese	x	x		online	
Portuguese, Brazilian	x			online	
Romanian	x			online	
Russian	x			online	
Serbian			male teen?	online	
Serbo-Croatian		x		online	
Sinhala	x			online	
Slovak	x			online	
Spanish	x			online	
Spanish, Latin American	x	x		online	
Swahili		x		online	

Swedish	x	x		online	
Tamil	x	x		online	
Thai	x	x		online	
Turkish	x	x		online	
Ukrainian	x			online	
Vietnamese	x	x		online	

## Windows Tablet

- Recommended browser: Microsoft Edge.
- Microsoft's online and offline voices are listed in the drop-down list of voices in AsTeRICS GRID. The quality of offline voices is generally lower than that of the online voices.
- While Microsoft often offers multiple versions for a language, differentiating by country, the voices do not always reflect the specific accent or dialect specific to that country!
- Microsoft does not explicitly label child voices. A few voices with female names sound very young. In the table, a cross was only marked if the voice distinctly sounded like a child. That said, several other "female" voices could also be used as child voices. A few of the voices are listed by Microsoft as "male" but don't sound deep. It might be a mistake by Microsoft or possibly a male teen before voice break.
- Voices installed through third-party software (e.g., CereProc) don't appear to work in the browser.

Language	Female voice	Male Voice	Child Voice	offline / online
Bulgarian (Bulgaria)		x		offline
Catalan (Spain)	x			offline
Chinese (Simplified PRC)	x			offline
Chinese (Traditional Hong Kong S.A.R.)	x	x		offline
Croatian (Croatia)		x		offline
Czech (Czech Republic)		x		offline
Danish (Denmark)	x			offline
Dutch (Belgium)		x		offline
English (Australia)	x	x		offline
English (Canada)	x	x		offline
English (Ireland)		x		offline
English (United Kingdom)	x	x		offline
English (United States)	x	x		offline
Finnish (Finland)	x			offline
French (Canada)	x	x		offline
French (France)	x	x		offline
French (Switzerland)		x		offline
German (Austria)		x		offline
German (Germany)	x	x		offline
German (Switzerland)		x		offline
Greek (Greece)		x		offline
Hebrew {Israel}		x		offline
Hindi (India)	x	x		offline
Hungarian (Hungary)		x		offline
Indonesian (Indonesia)		x		offline
Italian (Italy)	x	x		offline
Japanese (Japan)	x	x		offline

Malay (Malaysia)		x		offline
Norwegian (Bokmål Norway)		x		offline
Polish (Poland)	x	x		offline
Portuguese (Brazil)	x	x		offline
Portuguese (Portugal)	x			offline
Romanian (Romania)		x		offline
Russian (Russia)		x		offline
Slovak (Slovakia)		x		offline
Slovenian (Slovenia)		x		offline
Spanish (Mexico)	x	x		offline
Spanish (Spain)	x	x		offline
Swedish		x		offline
Thai (Thailand)		x		offline
Turkish (Turkey)		x		offline
Vietnamese (Vietnam)		x		offline
Afrikaans (South Africa)	x	x		online
Albanian (Albania)	x	x		online
Amharic (Ethiopia)	x	x		online
Arabic (Algeria)	x	x		online
Arabic (Bahrain)	x	x		online
Arabic (Egypt)	x	x		online
Arabic (Iraq)	x	x		online
Arabic (Jordan)	x	x		online
Arabic (Kuwait)	x	x		online
Arabic (Lebanon)	x	x		online
Arabic (Libya)	x	x		online
Arabic (Morocco)	x	x		online
Arabic (Oman)	x	x		online
Arabic (Qatar)	x	x		online
Arabic (Saudi Arabia)	x	x		online
Arabic (Syria)	x	x		online
Arabic (Tunisia)	x	x		online
Arabic (United Arab Emirates)	x	x		online
Arabic (Yemen)	x	x		online
Arabic, online	x	x		online
Armenian Male, online		x		online
Azerbaijani (Azerbaijan)	x	x		online
Bangla (Bangladesh)	x	x		online
Bangla (India)		x		online
Bangla India, online	x	x		online
Bengali (India)	x			online
Bosnian (Bosnia)	x	x		online
Bosnian, online		x		online

Bulgarian (Bulgaria)	x	x		online
Burmese (Myanmar)	x	x		online
Catalan (Spain)	x	x		online
Central Khmer (Cambodia)	x	x		online
Chinese (Cantonese Traditional)	x			online
Chinese (Hong Kong)	x			online
Chinese (Hong Kong)	x	x		online
Chinese (Mainland)	x	x	x	online
Chinese (Northeastern Mandarin)	x			online
Chinese (Taiwan)	x	x		online
Chinese (Taiwanese Mandarin)	x			online
Chinese (Zhongyuan Mandarin Shaanxi)	x			online
Chinese Male, online	x			online
Croatian (Croatia)	x	x		online
Czech (Czech Republic)	x	x		online
Danish (Denmark)	x	x		online
Dutch (Belgium)	x	x		online
Dutch (Netherlands)	x	x		online
English (Australia)	x	x		online
English (Canada)	x	x		online
English (Hongkong)	x	x		online
English (India)	x			online
English (Ireland)	x	x		online
English (Kenya)		x		online
English (New Zealand)	x			online
English (Nigeria)	x	x		online
English (Nigeria)		x		online
English (Philippines)	x	x		online
English (Singapore)	x	x		online
English (South Africa)	x	x		online
English (Tanzania)	x	x		online
English (United Kingdom)	x	x	x	online
English (United States)	x	x	x	online
Estonian (Estonia)	x	x		online
Filipino (Philippines)	x	x		online
Finnish (Finland)	x	x		online
French (Belgium)	x	x		online
French (Canada)	x	x		online
French (France)	x	x	x	online
French (Switzerland)	x	x		online
Galician (Spain)	x	x		online
Georgian (Georgia)	x	x		online



German (Austria)	x	x		online
German (Germany)	x	x		online
German (Switzerland)	x	x		online
Greek (Greece)	x	x		online
Gujarati (India)	x	x		online
Hebrew (Israel)	x	x		online
Hindi (India)	x	x		online
Hungarian (Hungary)	x	x		online
Icelandic (Iceland)	x	x		online
Indonesian (Indonesia)	x	x		online
Inuktitut (Canada)	x	x		online
Irish (Ireland)	x	x		online
Italian (Italy)	x	x		online
Japanese (Japan)	x	x		online
Javanese (Indonesia)	x	x		online
Kannada (India)	x	x		online
Kazakh (Kazakhstan)	x	x		online
Korean (Korea)	x	x		online
Lao (Laos)	x	x		online
Latin, online		x		online
Latvian (Latvia)	x	x		online
Lithuanian (Lithuania)	x	x		online
Macedonian (Republic of North Macedonia)	x	x		online
Malay (Malaysia)	x	x		online
Malayalam (India)	x	x		online
Maltese (Malta)	x	x		online
Marathi (India)	x	x		online
Moldavian, online	x			online
Mongolian (Mongolia)	x	x		online
Montenegrin, online	x			online
Nepali (Nepal)	x	x		online
Norwegian (Bokmål Norway)	x	x		online
Norwegian Male, online	x	x		online
Pashto (Afghanistan)	x	x		online
Persian (Iran)	x	x		online
Polish (Poland)	x	x		online
Portuguese (Brazil)	x	x		online
Portuguese (Portugal)	x	x		online
Romanian (Romania)	x	x		online
Russian (Russia)	x	x		online
Serbian (Serbia)	x	x		online
Sinhala (Sri Lanka)	x	x		online
Slovak (Slovakia)	x	x		online

Slovenian (Slovenia)	x	x		online
Somali (Somalia}	x	x		online
Spanish (Argentina)	x	x		online
Spanish (Bolivia)	x	x		online
Spanish (Chile)	x	x		online
Spanish (Colombia)	x	x		online
Spanish (Costa Rica)	x	x		online
Spanish (Cuba)	x	x		online
Spanish (Dominican Republic)	x	x		online
Spanish (Ecuador)	x	x		online
Spanish (El Salvador)	x	x		online
Spanish (Equatorial Guinea)	x	x		online
Spanish (Guatemala)	x	x		online
Spanish (Honduras)	x	x		online
Spanish (Mexico)	x	x		online
Spanish (Nicaragua)	x	x		online
Spanish (Panama)	x	x		online
Spanish (Paraguay)	x	x		online
Spanish (Peru)	x	x		online
Spanish (Puerto Rico)	x	x		online
Spanish (Spain)	x	x		online
Spanish (United States)	x	x		online
Spanish (Uruguay)	x	x		online
Spanish (Venezuela)	x	x		online
Spanish Latin American Female, online	x	x		online
Sundanese (Indonesia)	x	x		online
Swahili (Kenya)	x	x		online
Swahili (Tanzania)	x	x		online
Swedish (Sweden)	x	x		online
Tamil (India)	x	x		online
Tamil (Malaysia)	x			online
Tamil (Singapore)	x	x		online
Tamil (Sri Lanka)	x	x		online
Telugu (India)	x	x		online
Thai (Thailand)	x	x		online
Turkish (Turkey)	x	x		online
Ukrainian (Ukraine)	x	x		online
Urdu (India)	x	x		online
Urdu (Pakistan)	x	x		online
Uzbek (Uzbekistan)	x	x		online
Vietnamese (Vietnam)	x	x		online
Welsh (United Kingdom)	x	x		online
Zulu (South Africa)	x	x		online

## Android Tablet

- Recommended browser: Firefox.
- Android's online and offline voices are listed in the drop-down list of voices in AsTeRICS GRID. No voices are listed as offline, although some voices do appear to work offline.
- Voices purchased separately become available in the dropdown and are available offline. Tested with: Acapella, Vocalizer, RHVoice, Google Voices, AhoTTS, and CereProc. It does not seem possible to use voices from different TTS companies simultaneously. Set which TTS company to use under 'General management > Text-to-speech > Preferred engine.'
- Recommendation: Install voices for offline use. Tables with the voices available from the tested vendors are listed below.

Language	Female Voice	Male Voice
Afrikaans	x	
Albanian		x
Arabic	x	x
Bengali Bangla Bangladesh		x
Bengali India	x	x
Bosnian		x
Catalan	x	
Chinese	x	x
Chinese Hong Kong	x	x
Chinese Taiwan	x	x
Croatian	x	
Czech	x	
Danish	x	
Dutch	x	x
English Australia	x	x
English UK	x	x
English US	x	x
Estonian		x
Filipino	x	
Finnish	x	
French	x	x
French, Canada	x	x
German	x	x
Greek	x	
Hindi	x	x

Hungarian	x	
Icelandic	x	
Indonesian	x	x
Italian	x	x
Japanese	x	x
Korean	x	x
Latvian		x
Moldovian	x	
Montenegrin	x	
Nepali	x	
Norwegian Bokmål	x	x
Polish	x	x
Portuguese	x	x
Portuguese Brazil	x	
Romanian	x	
Russian	x	
Serbian		x
Serbo-Croatian	x	
Sinhalese	x	
Slovak	x	
Spanish	x	
Spanish LatAm	x	x
Swahili		x
Swedish	x	x
Tamil	x	x
Thai	x	x
Turkish		x
Ukrainian	x	
Vietnamese	x	x
Welsh	x	

Armenian and Macedonian are listed too but don't seem to function.

Several languages are labelled as male but sound female, and vice versa. The table reflects how the voices sound rather than their labels.

### **Acapela**

- Voices are purchased through the acapela app.

- Purchased voices are displayed in the drop-down, as e.g. “lang.spa: Spanish (ESP,lnes), online”. You need to select “show all voices”. Despite being labelled as “online” they are available offline too.
- All downloaded languages can be selected. If a voice is set as system voice it will appear as “Default device voice, offline” in the drop-down.

Language	Female Voice	Male Voice	Child Voice
Arabic	x	x	
Catalan	x		
Czech	x		
Danish	x		
Danish		x	
Dutch	x	x	boy, girl
Dutch, Belgium	x	x	
English, Australia	x	x	boy, girl
English, Canada	x		boy, girl
English, India	x		
English, Northern UK			child
English, UK	x	x	boy, girl
English, US	x	x	boy, girl teenage boy, teenage girl
Finnish	x		
French	x	x	boy, girl
French, Canada	x	x	boy, girl
German	x	x	boy, girl
Greek		x	
Hindi India	x		
Italian	x	x	boy, girl
Japanese	x		
Korean	x		
Mandarin Chinese	x		
Norwegian	x	x	boy, girl
Polish	x	x	
Portuguese	x		
Portuguese, Brazil	x	x	
Russian	x		
Scottish	x		boy, girl
Spanish	x	x	
Spanish, US	x	x	boy, girl

Swedish	x	x	boy, girl
Swedish, Finland		x	
Swedish, Göteborg		x	
Swedish, Scania	x		
Turkish	x		

Faroese is listed but doesn't seem to function.

### Vocalizer Voices

- Voices are purchased through the Vocalizer Voices app.
- Purchased voices are displayed in the drop-down as e.g. "lang.ron: Romanian (ROU), online". You need to select "show all voices". Despite being labelled as "online" they are available offline too.
- All downloaded languages can be selected.

Language	Female Voice	Male Voice	Child Voice
Arabic	x	x	
Basque	x		
Bengali, India			boy
Bhojpuri	x		
Bulgarian	x		
Cantonese (Hong Kong)	x	x	
Catalan	x	x	
Chinese, Mandarin	x	x	girl
Croatian	x		
Czech	x	x	
Danish	x	x	
Dutch	x	x	
English, Australia	x	x	
English, India	x	x	
English, Ireland	x		
English, Scotland	x		
English, South Africa	x		
English, UK	x	x	
English, USA	x	x	girl
Farsi		x	
Finnish	x	x	
Flemish	x		

French	x	x	
French, Belgium	x		
French, Canada	x	x	
Galician	x		
German	x	x	
Greek	x	x	
Hebrew	x		
Hindi	x	x	
Hungarian	x		
Indonesian	x		
Italian	x	x	
Japanese	x	x	
Kannada			girl
Korean	x	x	
Malay	x		
Mandarin Taiwanese	x		
Marathi			girl
Norwegian	x	x	
Polish	x	x	
Portuguese	x	x	
Portuguese, Brazil	x	x	
Romanian	x		
Russian	x	x	
Sichuanese	x		
Slowak	x		
Spanish	x	x	
Spanish Argentina	x	x	
Spanish Mexico	x	x	
Spanish, Chile	x		
Spanish, Columbia	x	x	
Swedish	x	x	
Tamil	x		
Telugu	x		
Thai	x		
Turkish	x	x	
Ukrainian	x		
Vietnamese	x		

Slowenian is listed but doesn't seem to function.

## RHVoice

- Voices are purchased through the RHVoice app.
- Purchased voices are displayed in the drop-down, as e.g. “Albanian (North Macedonia), online”. Despite being labelled as “online” they are available offline too.
- All downloaded languages can be selected.

Language	Female Voice	Male Voice
Albanian, North Macedonia	x	
Croatian	x	
Czech		x
Englisch, UK		x
English, India		x
English, Russia	x	x
English, USA	x	x
Georgian	x	
Kirgizh	x	x
Macedonian	x	x
Nepalese	x	
Polish	x	x
Portuguese, Brazil	x	x
Russian	x	x
Serbian	x	
Slovak	x	x
Spanish		x
Spanish, Mexico	x	
Tatar		x
Tswana	x	
Turkmenish	x	
Ukrainian	x	x
Uzbek	x	x
Vietnamese	x	

## Speech Recognition and Synthesis from google

- No purchase. Voices must be downloaded in the settings (Cog wheel > Install voice data). If multiple voices are available for the same language, changing it in the settings doesn't seem to impact the voice that is used in the browser. The table shows the voices that are automatically selected.



- The voices are displayed in the drop-down, as e.g. “Khmer (Cambodia), offline”.
- All downloaded languages can be selected.
- There is sometimes a lag when switching languages
- While Android often offers multiple versions for one language, differentiating by country, the voice is most of the time just the same voice.

Language	Female Voice	Male Voice
Albanian (Albania)		x
Arabic		x
Bagla (Bangladesh)	x	
Bangla (India)	x	
Bosnian (Bosnia & Herzegovina)		x
Bulgarian (Bulgaria)	x	
Cantonese (Hong Kong)	x	
Catalan (Spain)	x	
Central Khmer (Cambodia)	x	
cmn (China)	x	
cmn (Taiwan)	x	
Croatian (Croatia)	x	
Danish (Denmark)	x	
Dutch (Belgium)	x	
Dutch (Netherlands)	x	
English (Australia)	x	
English (India)	x	
English (Nigeria)	x	
English (United Kingdom)	x	
English (United States)	x	
Estonian (Estonia)		x
Filipino (Philippines)	x	
Finnish (Finland)	x	
French (Canada)	x	
French (France)	x	
German (Germany)	x	
Greek (Greece)	x	
Gujarati (India)	x	
Hebrew (Israel)	x	
Hindi (India)	x	
Hungarian (Hungary)	x	
Icelandic (Iceland)	x	

Indonesian (Indonesia)	x	
Italian (Italy)	x	
Japanese (Japan)	x	
Javanese (Indonesia)	x	
Kannada (India)	x	
Korean (South Korea)	x	
Latvian (Latvia)		x
Lithuanian (Lithuania)		x
Malay (Malaysia)	x	
Malayalam (India)	x	
Marathi (India)	x	
Nepali (Nepal)	x	
Norwegian Bokmål (Norway)	x	
Polish (Poland)	x	
Portuguese (Brazil)	x	
Portuguese (Portugal)	x	
Punjabi (India)	x	
Romanian (Romania)	x	
Russian (Russia)	x	
Sindhi (India)		
Sinhala (Sri Lanka)	x	
Slovak (Slovakia)	x	
Slovenian (Slovenia)	x	
Spanish (Spain)	x	
Spanish (United States)	x	
Sundanese (Indonesia)	x	
Swahili (Kenya)		x
Swedish (Sweden)	x	
Tamil (India)	x	
Telugu (India)	x	
Thai (Thailand)	x	
Turkish (Turkey)	x	
Ukrainian (Ukraine)	x	
Urdu (Pakistan)		x
Vietnamese (Vietnam)	x	
Welsh (United Kingdom)	x	

The following languages don't seem to function: Assamese (India), Bodo (India), Dogri (India), Kashmiri (India), Konkani (India), Maithili (India), Manipuri (India), Odia (India), Sanskrit (India), Santali (India); Czech (Czechia)

**AhoTTS**

Offers voices in Spanish and Basque

**CereProc**

- Only the system voice seems available as default.
- It doesn't work with all purchased voices.

→ Not recommended

## Voices Per System

*Only languages that available offline are listed. The quality varies.*

<b>iOS, Windows, Android</b>	<b>iOS, Android</b>	<b>Android</b>	<b>iOS</b>
Cantonese	Arabic	Albanian	Kazakh
Catalan	Bangla	Bengali	Persian
Chinese (Mandarin)	Basque	Bosnian	Shanghainese
Croatian	Bhojpuri	Central Khmer	
Czech	Bulgarian	Estonian	
Danish	Galician	Farsi	
Dutch	Kannada	Filipino	
Englisch	Korean	Flemish	
Finnish	Lithuanian	Georgian	
French	Marathi	Gujarati	
German	Taiwan	Icelandic	
Greek	Tamil	Kirgizh	
Hebrew	Telugu	Latvian	
Hindi	Ukrainian	Macedonian	
Hungarian		Malayalam	
Indonesian		Nepali	
Italian		Punjabi	
Japanese		Scottish	
Malay		Serbian	
Norwegian, Bokmål		Sichuanese	
Polish		Sindhi	
Portuguese		Sinhala	
Romanian		Slovak	
Russian		Sundanese	
Slovak		Swahili	
Slovenian		Taiwanese	
Spanish		Tatar	
Swedish		Tswana	
Thai		Turkmenish	
Turkish		Urdu	
Vietnamese		Uzbek	
		Welsh	